



UD Sleeve Sampler



Volume	Length mm (A)
0.1cc	63
0.2cc	61
0.25cc	59
0.3cc	57.5
0.4cc	54
0.5cc	50
0.6cc	46.5
0.7cc	43
0.8cc	39
0.9cc	35
1cc	31
1.1cc	27
1.2cc	23

Volume	Length mm (B)
1.3cc	42
1.4cc	40
1.5cc	38
1.6cc	36
1.7cc	34
1.8cc	32
1.9cc	30
2.0cc	28
2.1cc	26
2.2cc	24
2.3cc	22
2.4cc	20
2.5cc	18