

Certificate of Conformity

Part No.	8063H-300 8063H-300S 8063H-500 8063H-500S 8063H-1000 8063H-1000S
Description	Powder Lance
Material	HDPE (High Density Polyethylene)
Date	13-Sep-10

Dear Sir/Madam,

We confirm that the goods are manufactured with HDPE that conforms to FDA CFR 177.1520.

The samplers are assembled and packed in a class 100,000 medical cleanroom.

Ultra-sonic welding is used adhere the components (no adhesives, screws etc. are used).

The HDPE resin used in the manufacture of the above products is BSE/TSE free.

For and on behalf of
Sampling Systems,



Jonathan Fenner